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**Robert Barrack, MD**  
St. Louis, Missouri

**Robert Bourne, MD**  
London Health Sciences Center  
London, Ontario, Canada

**Jonathan Garino, MD**  
University of Pennsylvania School of Medicine  
Philadelphia, Pennsylvania

**Wayne M. Goldstein, MD**  
Clinical Professor of Orthopaedics  
University of Illinois at Chicago  
Illinois Bone and Joint Institute  
Chicago, Illinois

**Richard Kyle, MD**  
Minneapolis, Minnesota

**Stephen J. McMahon MB BS,  
FRACS(Orth), FA(Orth)A**  
Senior Lecturer Monash University  
Malabar Orthopaedic Clinic  
Melbourne, Australia

**John L. Masonis, MD**  
OrthoCarolina  
Hip & Knee Center  
Charlotte, North Carolina

**Henrik Malchau, MD**  
Associate Professor Harvard Medical School  
Codirector The Harris Orthopaedic  
Biomechanics and Biomaterials Laboratory  
Massachusetts General Hospital  
Boston, Massachusetts

**Michael Ries, MD**  
University of California  
San Francisco, California

**Cecil Rorabeck, MD**  
Professor of Orthopaedic Surgery  
University of Western Ontario  
London, Ontario, Canada

**Van Paul Stamos, MD**  
Illinois Bone and Joint Institute  
Glenview, Illinois  
Clinical Instructor of Orthopaedic Surgery  
Northwestern University Medical School  
Chicago, Illinois

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## Referencias

- 1 Biolox-ceramics for his arthroplasty. CeramTec AG, MT 060003: GB.5.000-0612. Germany.
- 2 CeramTec Technical Monograph. Current perspective on the use of ceramics in total hip arthroplasty. CeramTec AG, 060003: GB.5.000-0612, 2007.
- 3 Internal Smith & Nephew testing on file.
- 4 Li MG, Zhou ZK, Wood DJ, et al. Low wear with high-cross linked polyethylene especially in combination with OXINIUM heads. A RSA evaluation. Poster no. 643 presented at: Orthopaedic Research Society Annual Meeting; March 19–22, 2006; Chicago, IL.
- 5 Clarke IC, Good P, Williams P, et al. Ultra-low wear rates for rigid-on-rigid bearings in total hip replacements. *Poc Inst Mech Eng. [IJ]* 2000;214(4):331–347.
- 6 Good V, Ries m, Barrack RL, Widding K, Hunter G, Heuer D. Reduced wear with oxidized zirconium femoral heads. *J Bone Joint Surg Am.* 2003;85-A Suppl 4:105–110.
- 7 Scott M, Morrison M, Mishra SR, Jani S. A method to quantify wear particle volume using atomic force microscopy. *ORS Transactions.* 2002;27:132.
- 8 Ries MD, Scott ML. Relationship between gravimetric wear and particle generation in hip simulators: conventional versus crosslinked polyethylene. Scientific exhibit at American Academy of Orthopaedic Surgeons; Feb 27–March 4, 2001; San Francisco, CA.
- 9 Good V, Widding K, Heuer D, Hunter G. Reduced wear using the ceramic surface on oxidized zirconium heads. In: Lazennec JY, Dietrich M, eds. *Bioceramics in Joint Arthroplasty.* Darmstadt, Germany: Steinkopff; 2004:93–98.
- 10 Berry DJ, von Knoch M, Schleck CD, Harmsen WS. Effect of femoral head diameter and operative approach on risk of dislocation after primary total hip arthroplasty. *J Bone Joint Surg Am.* 2005 Nov;87(11):2456–2463.
- 11 Barrack RL, Butler RA, Laster DR, Andrews P. Stem design and dislocation after revision total hip arthroplasty: clinical results and computer modeling. *J Arthroplasty.* 2001 Dec;16(8 Suppl 1):8–12.
- 12 Barrack RL. Dislocation after total hip arthroplasty: implant design and orientation. *J Am Acad Orthop Surg.* 2003 Mar-Apr;11(2):89–99.
- 13 Barrack RL, Lavernia C, Ries M, Thornberry R, Tozakoglou E. Virtual reality computer animation of the effect of component position and design on stability after total hip arthroplasty. *Orthop Clin North Am.* 2001 Oct;32(4):569–577, vii.
- 14 Tradonsky S, Postak P, Frimson A, Greenwald A. Performance characteristics of two piece acetabular cups. Cleveland, OH: The Orthopaedic Research Laboratory, Mt. Sinai Medical Center. 1992.
- 15 FDA guidance document for testing acetabular cup prosthesis. Orthopedic Services Branch, US Food and Drug Administration. May 1995.
- 16 Bourne R. Randomized controlled trial to compare acetabular component fixation of two porous ingrowth surfaces using RSA analysis. London, Ontario, Canada: London Health Science Center. 2007. Internal report on file at Smith & Nephew, Memphis, TN.
- 17 Heiner AD, Brown TD. Frictional coefficients of a new bone ingrowth structure. Poster no. 1623 presented at: Orthopaedic Research Society Annual Meeting; Feb 11–14, 2007; San Diego, CA.
- 18 Padgett DE, Miller AN, Su EP, Bostrom MPG, Nestor BJ. Ceramic liner malseating in total hip arthroplasty. Poster PO97 at American Academy of Orthopaedic Surgeons; Feb 14–18, 2007; San Diego, CA.
- 19 Langdown AJ, Pickard RJ, Hobbs CM, Clarke HJ, Dalton DJ, Grover ML. Incomplete seating of the liner with the Trident acetabular system: a cause for concern? *J Bone Joint Surg Br.* 2007 Mar;89(3):291–295.

**Divisi3n ASD**  
**Smith & Nephew, S.A.U.**  
Fructu3s Gelabert, 2-4  
08970 Sant Joan Desp3  
Barcelona (Espana)

T +34 93 373 73 01  
F +34 93 373 74 53  
[www.smith-nephew.com](http://www.smith-nephew.com)